



PurMa Biologics, LLC

## PurMa™ Mycoplasma Detection Kit

Cat #: P5M011404

### User Manual

#### The Kit Contents (50 Reactions) Cat # P5M011404-50:

- 1) PurMa™ PCR 2X Mix (containing green dye and includes all the primers) (700 µl)
- 2) Reaction Buffer (2500µl)
- 3) Positive control. The lysed mycoplasma which does not carry any live mycoplasma (50ul)
- 4) Negative Control. Lysed COS1 cells serves as the negative control (mycoplasma free), ready to proceed to the PCR and (50µl).
- 5) PurMa™ Ultra-Pure Nuclease Free Water (Cat# P3W013104) (1ml)
- 6)




Figure1. PurMa™ Mycoplasma Detection Kit produces a 230bp fragment.

Lanes:

1. PurMa™ Negative Control.
2. PurMa™ Positive Control.
3. Contaminated medium.
4. Contaminated FBS.
5. Contaminated Cell line.

## **Protocol**

### **Part 1: Tissue/Cell Homogenization**

#### Preparing Samples

- 1) Take 1000 µl of cell culture supernatant including the cells. This is because the aim is to detect all the hidden/ intracellular mycoplasma as well. place them in a 1.5ml tube.
- 2) Centrifuge at 13000 RPM for 10 min, discard the sup and add 50µl of PurMa™ Reaction Buffer.
- 3) Heat the sample for exactly 5min at 95°C (to desaturate the nuclear cell wall and release the mycoplasma from heavy membrane).

#### **Important Note:**

**Do not incubate for more than five min as there will be false positive result otherwise.**

- 4) Vortex for 20 seconds
- 5) Transfer immediately on the ice and incubate for 5 min.
- 6) During the incubation on ice, vortex them two more and put back on the ice immediately.
- 7) Spin down briefly
- 8) Take 1 µl of each sample for the PCR reaction.

### **Part 2: PCR Reaction**

PCR reaction Ingredients	µl
PurMa™ 2X PCR MIX	12.5
Sample from Part 1	1
PurMa™ Ultra-Pure Water	
Total	25

We strongly recommend preparing master mix and do not prepare tubes individually.

### **Part 3: Thermal Profile**

Perform the PCR using the following program:

Step	Temperature (°C)	Time
Initial denaturation	94	5 min
34cycles	94	30 S
	55	45 S
	68	1 min
Final extension	68	10 min
	4	Over night

### **Part 4: Running the result:**

Prepare a 2 % agarose and the product will show at 230 bp.

#### **Important Note:**

Due to mixing several primers, in some instances you might see a primer dimer band at 50 bp, please disregard that band as that has nothing to do with the presence/ absence of mycoplasma in your sample.